



TECHNICAL BULLETIN

Bureau of Clinical Laboratories, Sharon P. Massingale, Ph.D., HCLD (ABB) Director

September 15, 2014

Volume V Issue 4

ATTENTION HEALTH CARE PROVIDERS



When collecting specimens for CT/GC/TV please be mindful when selecting the specimen source. There are 3 collection kits available:

- **Aptima Unisex Swab Collection Kit:** Order # **CF02** Specimen source is Endocervical, Male Urethral, Pharyngeal (Oral) and Rectal.
- **Aptima Urine Collection Kit:** Order# **CM01** Male and Female Urine specimens
- **Aptima Vaginal Swab Collection Kit:** Order # **CF04** Vaginal swabs

Make sure the proper collection kit is used for the specimen collected. These Aptima Kits are available to order on ProcureIT by the order #s given above.

Please contact Traci Dailey at 334-260-3432 or Carol Dorsey at 251-344-6049 if you have questions regarding this bulletin.

*Bureau of Clinical Laboratories
Quality Management
8140 AUM Drive
Montgomery, Alabama 36117*

